

## BIOFREEZER SYSTEM OPTIONS

### Cryobath Options

#### LNS

1.5 liter insulated and covered metal canister LN2 cryobath for freezing cycles of up to 90 minutes.

PRODUCT CODE: LNS



#### LNL

3.8 liter insulated and covered metal canister LN2 cryobath for freezing cycles of up to 3.5 hours.

PRODUCT CODE: LNL



#### DX2

2.3 liter EPS foam cryobath with PETG clear lid for increased LN2 economy. For freezing cycles up to 2.25 hours, eliminating the need for refills during a freeze. Vent holes allow easy access for seeding instruments.

PRODUCT CODE: DX2

