

Labeling Systems

Brady® TLS2200 Labeling System

The TLS2200 is an efficient, cost-effective, and time-saving solution to the challenge of labeling biological specimens for cryostorage. This unique handheld instrument generates labels with specially formulated adhesive for straw ID rods, cryovials, test tubes, and other labware. The user has control of the type size, number of lines, orientation, and virtually every parameter of the print format. Bar coding is also supported by the TLS2200.

The TLS2200 comes standard with attaché-style carry case, rechargeable battery and AC adapter (specify USA or international version), sample labels, protective cover, cleaning kit, and comprehensive instructions.

PRODUCT CODE: TLS

