

## Microscopes

### [Wesco 701 Zoom Stereo Microscope](#)

The 701 Zoom Stereo Microscope is designed for superior optical and mechanical performance. With its wide 7:1 zoom ratio and choice of accessories, the 701 is an extremely versatile instrument. Many configurations are available, making the 701 the clear choice for industrial or bio-medical use.

#### 701 Specifications

Working distance 100mm

Interpupillary distance 55-75mm

Diopter adjustment  $\pm 5$

Focus range 135mm

Observation body 45° inclined, 360° rotatable

Stage 95mm stage plate

PRODUCT CODE: W701



### [Wesco CXRIII Compound Microscope](#)

WESCO CXRIII is a laboratory microscope system for high-end application, with a modern design incorporating the latest optical and mechanical advancements. Extra clarity and contrast is provided through high performance infinity corrected IP optical system. WESCO CXRIII offers great performance features standard, making it a clear choice for advanced laboratory and research applications.

PRODUCT CODE: WCXR

